


4/8/2009

MITK-IGT

Introduction to the image guided therapy component
of the Medical Imaging Interaction Toolkit

Intra-interventual visualization of instruments in relation to patient using preoperative or intra operative images



MITK

Application Framework

GUI

Rendering

Interaction

Data
Management


Algorithms

PACS

IGT


Base Libraries




- NDI Polaris
- NDI Aurora
- Claron Technologies MicronTracker (Windows only)
- Ascension Microbird (Windows only)
- Virtual tracking device




- Supports Windows, Linux (and MacOS) platforms
- Seamless integration with MITK
- Manages communication with tracking devices
- Provides Hardware abstraction of tracking devices
- Processing of tracking data with Pipeline architecture
- Ready to use modules for:
 - Recording and Replaying of tracking data
 - Coordinate transformations, including Landmark based transformation
 - Synchronization of multiple tracking devices
 - visualization of instrument movement
 - Visualization of camera movement
 - Visualization of movement trajectories
- Easily extensible


Example MITK-IGT Filter Pipeline








IGT tracking device classes




IGT utility classes



- Aim of the example:
 - Graphical application
 - read tracking data from tracking device
 - visualize movement of tracked object

- Tutorial is based on (old) Qt3 MainApp
- Create Functionality:



- Move to /mitk/QFunctionalities/

- Edit /mitk/QFunctionalities/QmitkBVMIGTTutorial/CMakeLists.txt

```
CREATE_QFUNCTIONALITY(QmitkBVMIGTTutorial)
INCLUDE_DIRECTORIES(${MITK_IGT_INCLUDE_DIRS})
LINK_LIB_TO_FUNCTIONALITY(QmitkBVMIGTTutorial mitkIGT)
```

- Activate Functionality using Cmake


- `BUILD_QmitkBVMIGTTutorial` must be set to ON
- `MITK_BUILD_IGT` must be set to ON

- Re-generate MITK-Project files using Cmake
- load MITK-Project in development environment

- Add members to class definition:

```
class QmitkIGTTutorial : public QmitkFunctionality
{
    [...]
protected slots:
    void OnSetupIGT();
    void OnTimer();
    void OnStop();
protected:
    mitk::TrackingDeviceSource::Pointer m_Source;
    mitk::NavigationDataVisualizationByBaseDataTransformFilter
        ::Pointer m_Visualizer;
    QTimer* m_Timer;
}
```

- Add Start and Stop buttons to user interface
- Connect buttons to methods



```
void QmitkBVMIGTTutorialTest::CreateConnections ()
{
    if (m_Controls)
    {
        connect ((QObject*) (m_Controls->m_StartButton),
                SIGNAL(clicked()), (QObject*) this, SLOT (OnSetupIGT ()));
        connect ((QObject*) (m_Controls->m_StopButton),
                SIGNAL(clicked()), (QObject*) this, SLOT (OnStop ()));
    }
}
```

- Implement Setup method: Setup tracking device

```
void QmitkIGTTutorial::OnSetupIGT()
{
    mitk::NDITrackingDevice::Pointer tracker =
        mitk::NDITrackingDevice::New();
    tracker->SetPortNumber(mitk::SerialCommunication::COM4);
    tracker->SetBaudRate(mitk::SerialCommunication::BaudRate115200);
    tracker->SetType(mitk::NDIPolaris);

    mitk::NDIPassiveTool::Pointer tool = mitk::NDIPassiveTool::New();
    tool->SetToolName("MyInstrument");
    tool->LoadSRROMFile("c:\\myinstrument.rom");
    tracker->Add6DTool(tool);
    ...
}
```


- Implement Setup method: Setup IGT filter pipeline

```
...
m_Source = mitk::TrackingDeviceSource::New();
m_Source->SetTrackingDevice(tracker);
m_Source->Connect();

m_Visualizer =
    mitk::NavigationDataVisualizationByBaseDataTransformFilter::New();
m_Visualizer->SetInput(m_Source->GetOutput());
mitk::Cone::Pointer instrument = mitk::Cone::New();
m_Visualizer->SetBaseData(m_Source->GetOutput(), instrument);

mitk::DataTreeNode::Pointer node = mitk::DataTreeNode::New();
node->SetData(instrument);
node->SetName("My tracked object");
node->SetColor(1.0, 0.0, 0.0);
mitk::DataStorage::GetInstance()->Add(node);
...
```

- Implement Setup method: start continuous update

...

```
m_Source->StartTracking();  
m_Timer = new QTimer(this);  
connect(m_Timer, SIGNAL(timeout()), this, SLOT(OnTimer()));  
m_Timer->start(100);  
}
```

- Implement Update method: execute continuous update

```
void QmitkIGTTutorial::OnTimer()  
{  
    m_Visualizer->Update();  
    mitk::RenderingManager::GetInstance()->RequestUpdateAll();  
}
```

- Implement Stop method: cleanup everything

```
void QmitkIGTTutorial::OnStop()  
{  
    m_Timer->stop();  
    disconnect(m_Timer, SIGNAL(timeout()), this, SLOT(OnTimer()));  
    m_Source->StopTracking();  
    m_Source->Disconnect();  
}
```

- Add include files to header:

- `#include "mitkTrackingDeviceSource.h"`
- `#include "mitkNavigationDataVisualizationByBaseDataTransformFilter.h"`
- `#include "qtimer.h"`

- Add include files to .cpp:

```
#include "mitkNDITrackingDevice.h"  
#include "mitkNDIPassiveTool.h"  
#include "mitkCone.h"
```

- Compile
- Execute
- Navigate!

- MITK-IGT will be available for download within the next days at:
`svn co http://svn.mitk.org/`
- A Software demo will be held at **Monday 16:40** in the Foyer
- More documentation is available:
 - BVM Proceedings Paper **1399**: *MITK-IGT: Eine Navigationskomponente für das Medical Imaging Interaction Toolkit*
 - MITK Documentation: <http://mitk.org/wiki/Documentation>
 - 2 Tutorials and an example application module available:
 - `/mitk/Core/IGT/IGTTutorial/IGTTutorialStep1.cpp`
 - `/mitk/Qfunctionalities/QmitkIGTTutorial`
 - `/mitk/Qfunctionalities/QmitkIGTExample`